

Available on FP10/GP10 - PLEASE COMPLETE ALL SECTIONS
TO BE COMPLETED BY THE HEALTHCARE PROFESSIONAL AND GIVEN TO YOUR PRESCRIBER

1 Patient Name

Doctor/Clinician Clinic

Pharmacist

Date

2 KITS:

Contain 2 x White Closed Toe Liners (10 mmHg each) and
1 x Outer Open Toe Stocking (30 mmHg) in Soft Beige

QUANTITY:

LINER PACKS:

Liner packs contains 3 white liners providing
10mmHg each

COLOUR:

Soft Beige White


QUANTITY:

3 MEASUREMENTS:

MEASUREMENT TIPS:

- Measurements should be taken with the person standing, feet flat on the floor. If your patient is unable to stand, you can take the measurements with the person seated, and feet flat on the floor.
- For foot measurement: stand on a plain sheet of paper, mark at the toe and at the heel; remove the foot and measure between the two points; this will provide a more accurate reading.
- When possible, measurements should be taken first thing in the morning.⁽¹⁾

1. NICE Clinical Knowledge Summary - Compression Stockings (2008, last update 2022)



All measurements are in centimetres

	B: Calf	C: Ankle	Size Required
Small	32 - 38	19 - 25.5	
Medium	34 - 40	21 - 27.5	
Large	36 - 42	23 - 29.5	
X Large	38 - 45	25 - 31.5	
XX Large	40.5 - 50	26.5 - 34	

	L: Leg	Size Required (please tick)
Short	34 - 38	
Regular	38 - 42	
Long	42 - 46	

COMMENTS:

NEED HELP? Contact Urgo Medical for CPD accredited training. **T: 0330 128 0898**