

Available on FP10/GP10 - PLEASE COMPLETE ALL SECTIONS
TO BE COMPLETED BY THE HEALTHCARE PROFESSIONAL AND GIVEN TO YOUR PRESCRIBER

1 Purchase Order No.
Patient Name
Altiven Registration No. (Repeat Orders)

DELIVER TO PHARMACY/HOSPITAL:
Name
Address
..... Postcode

Telephone No.

MEASUREMENTS COMPLETED BY:
Health Centre

Contact No.

2 INVOICE TO:

Hospital/Clinic
 Buying Group
 Wholesaler
 Pharmacy Chain

Name of buying group

Account No.

3 COMPRESSION CLASS: Low (14-18 mmHg)
 1 (18-21 mmHg) 2 (23-32 mmHg) 3 (34-46 mmHg)

THIGH LENGTH STYLE OPTIONS: Standard Hold Ups - Plain
 Standard with waist attachment Hold Ups - Lace

OTHER STYLES: Below Knee Tights

FOOT: Closed Toe Open Toe If you have chosen thigh length stockings, choose what styles you would like, at no extra cost.

4 QUANTITY: **COLOUR:**


STOCKINGS (SINGLES): **TIGHTS (PAIRS):** Beige
Left Leg Right Leg Black

5 MEASUREMENTS:

MEASUREMENT TIPS:

- Measurements should be taken with the person standing, feet flat on the floor. If your patient is unable to stand, you can take the measurements with the person seated, and feet flat on the floor.
- For foot measurement: stand on a plain sheet of paper, mark at the toe and at the heel; remove the foot and measure between the two points; this will provide a more accurate reading.
- When possible, measurements should be taken first thing in the morning.⁽¹⁾
- For thigh length, take the measurement while the person is standing.

1. NICE Clinical Knowledge Summary - Compression Stockings (2008, last update 2022)



Please fill in all the appropriate measurements in centimetres

CIRCUMFERENCE MEASUREMENTS	Left Leg	Right Leg
CT (waist - for tights only)		
CH (hips - for tights only)		
CG (thigh)		
CD (under knee)		
CC (calf)		
CB (ankle)		
CA (measure at base of little toe)		

LENGTH MEASUREMENTS	Left Leg	Right Leg
LT (for tights only)		
LG		
LD		
LA (for open toe only)		
LZ (for closed toe only)		

COMMENTS:

72 hours
despatch

SEND US YOUR FORM:
F: 01509 320 139 E: ordersin@uk.urgo.com
T: 0330 128 0898

WANT TO ORDER ONLINE?
www.altimed.co.uk

Doc Ref: 191346
Version: 02 Rev: 2024-02

NEED HELP? Contact Urgo Medical
for CPD accredited training.