



# COMPRESSION GARMENT RANGE



Accredits Urgo Medical research in  
wound healing and healthy skin

Altipress® Altiform® Altiven®


MADE IN  
BRITAIN 

ABOUT



## FULL COMPRESSION PORTFOLIO

Urgo Medical, a market leader in advanced woundcare and compression, presents a full compression portfolio including dual-layer bandages (**UrgoKTwo**), four-layer bandages (**KFour**) and hosiery, for all your patient's needs, from treatment to prevention.



## SUPPORT & TRAINING AT EVERY STEP

Urgo Medical regional teams are available to train, educate, hold study events and support the implementation of pathways.

- Hosiery study days, including CPD modules, practical sessions and product demonstrations
- Online resources for pharmacists
- Support designing, printing and delivering pathways.



## MADE IN THE UK

All our hosiery is manufactured in our factory in Shepshed, Loughborough. We offer site visits to demonstrate our commitment to the highest quality compression solutions.



## FAST DELIVERY

We manufacture locally, resulting in market-leading delivery time.

Once your Made-to-Measure order is received, it is manufactured within 72 hours.

# Altipress<sup>®</sup>

LEG ULCER KITS AND LINER PACKS

# Altiform<sup>®</sup>

BRITISH STANDARD COMPRESSION HOSIERY

# Altiven<sup>®</sup>

EUROPEAN STANDARD COMPRESSION HOSIERY

# Sockaid<sup>™</sup>

HOSIERY APPLICATOR



The Altipress 40 Leg Ulcer Hosiery kit offers a practical and discreet alternative to bulky bandage systems. It is ideal for a patient with the will and ability to self-care.

Consists of 2 stockings that are layered on top of each other:

- Layer 1: Liner - Applies 10 mmHg at the ankle
- Layer 2: British Standard Class III stocking - Applies 30 mmHg at the ankle

Together, they provide 40 mmHg at the ankle, which is required for leg ulcer healing. Altipress Liners are also available separately.

### BENEFITS OF ALTIPRESS:

- ✓ 15 Ready-to-wear sizes - options for all your patient's needs
- ✓ Easy to apply with a unique elasticated toe and heel design
- ✓ Ideal for self-caring patients
- ✓ British Skin Foundation Accredited - Demonstrating our commitment to quality products that are not harmful to the skin
- ✓ Latex Free

### INDICATIONS

- Treatment and management of venous leg ulcers
- Prevention of venous leg ulcer recurrence
- Patients not suitable for or tolerant of compression bandages

### CONTRAINDICATIONS

- Peripheral Vascular Disease (PVD), confirmed with an ABPI (Doppler) reading of less than 0.8
- Acute or unstable Heart Failure
- Ischaemic Heart Disease (IHD), also called Coronary Artery Disease
- Untreated Cellulitis or Phlebitis of the Lower Limb

### Proceed with caution in those patients who have:

- Diabetes
- Rheumatoid Arthritis
- Heavily exuding or large wounds
- Irregular shaped limb

These patients may require referral to a specialist for advice and additional assessment. Close supervision may be necessary in line with local policies/pathways.

### AVAILABLE IN:

#### STANDARD SIZES

Sizes: S, M, L, XL, XXL.

Lengths: Short, Regular, Long.

Kits contain 2 liners & 1 outer stocking

- Liner ○
- Stocking ●

#### MADE-TO-MEASURE

Tailor-made to patient's individual dimensions.

Kits contain 2 liners & 1 outer stocking

- Liner ○
- Stocking ●

Altipress Liners are available separately. Altipress Liners provide light compression (10 mmHg at the ankle) for maintenance therapy, palliative care and intolerance and/or non-concordance to higher levels of compression. They can be layered to provide higher levels of compression (e.g. 2 liners provide 20 mmHg at the ankle) or worn under other compression stockings.

## CONTRAINDICATIONS

If any existing arterial conditions exist, please seek advice from a healthcare professional before use.

## AVAILABLE IN:

### STANDARD SIZES

Sizes: S, M, L, XL, XXL

Lengths: Short, Regular, Long

Standard size liner packs contain 3 liners in either:

- White ○
- Soft Beige ●

## MADE-TO-MEASURE

Tailor-made to patient's individual dimensions.

Made-to-measure liner packs contain 2 liners in:


- White ○

## MEASURING AND FITTING GUIDE

	B: Calf	C: Ankle
Small	32 - 38	19 - 25.5
Medium	34 - 40	21 - 27.5
Large	36 - 42	23 - 29.5
X Large	38 - 45	25 - 31.5
XX Large	40.5 - 50	26.5 - 34

	Leg Length
Short	34 - 38
Regular	38 - 42
Long	42 - 46

Measurements in cms



To order our Made-to-Measure hosiery, please request an order pad or order online at [www.altimed.co.uk](http://www.altimed.co.uk)

**Altiform compression hosiery is specially designed to improve blood circulation in the legs. It provides graduated compression that gently squeezes the leg, pushing the blood back towards the heart and reducing fluid build-up at the ankle.**

**Altiform is recommended to prevent venous disease progression.**



### BENEFITS:

- ✓ Discreet, comfortable and easy to wear
- ✓ Exceptional stretch, making it easy to slip over the heel
- ✓ Unique elasticated toe & heel for extra comfort
- ✓ Latex Free
- ✓ British Skin Foundation Accredited - Demonstrating our commitment to quality products that are not harmful to the skin
- ✓ Made in the UK.

### INDICATIONS

Available in 3 Classes of Graduated Compression:

#### Class I – Light support

(14 – 17 mmHg compression at the ankle)

- Superficial or early varices
- Varicosis during pregnancy
- Swollen or aching legs & ankles
- Suitable for long distance travel

#### Class II – Medium support

(18 – 24 mmHg compression at the ankle)

- Varices of medium severity
- Venous leg ulcer treatment and prevention of recurrence
- Mild oedema
- Varicosis during pregnancy
- Suitable for long distance travel

#### Class III – Strong support

(25 – 35 mmHg compression at the ankle)

- Gross varices
- Post thrombotic venous insufficiency
- Gross oedema
- Venous leg ulcer treatment
- Prevention of recurrence

### CONTRAINDICATIONS


- Peripheral Vascular Disease (PVD), confirmed with an ABPI (Doppler) reading of less than 0.8
- Acute or unstable Heart Failure
- Ischaemic Heart Disease (IHD), also called Coronary Artery Disease
- Untreated Cellulitis or Phlebitis of the Lower Limb

#### Proceed with caution in those patients who have:

- Diabetes
- Rheumatoid Arthritis

These patients may require referral to a specialist for advice and additional assessment. Close supervision may be necessary in line with local policies/pathways.

## STYLES AVAILABLE:



BELOW KNEE




THIGH LENGTH



HOLD UPS PLAIN



HOLD UPS LACE



## AVAILABLE IN STANDARD SIZES: SMALL, MEDIUM, LARGE, X-LARGE

	A: Mid-thigh	B: Calf	C: Ankle	D: Foot
Small	41 - 54	32 - 38	19 - 25.5	20 - 26
Medium	44 - 58	34 - 40.5	21 - 27.5	22 - 28
Large	47 - 61	36 - 43.5	23 - 29.5	24 - 31
XLarge	50 - 64	38 - 47	25 - 32	26 - 33

Measurements in cms

	Below Knee		Thigh	
	Open Toe	Closed Toe	Open Toe	Closed Toe
Class I	●	● ●	●	● ●
Class II	●	● ●	●	● ●
Class III	●	●	●	●

Soft Beige ● Black ●

## AND MADE-TO-MEASURE FOR A TAILORED FIT

Well-fitting compression stockings make all the difference to a patient's well-being and can improve concordance. Altiform Made-to-Measure stockings are made individually for each leg - giving the patient a precise, comfortable fit. Right and left leg measurements are often different - so different that sometimes two pairs of stock are needed to account for these differences. For this reason, Made-to-Measure stockings are a cost-effective option.

## BENEFITS OF ALTIFORM MADE TO MEASURE HOSIERY:

- ✓ **Manufactured within 72 hours - ensuring fast delivery**
- ✓ Tailor-made to provide a unique and comfortable fit, to help improve concordance
- ✓ Easy to apply with a unique elasticated toe and heel design
- ✓ Unique patient identification number sewn into the sole for easy re-ordering
- ✓ Latex Free



To order our Made-to-Measure hosiery, please request an order pad or order online at [www.altimed.co.uk](http://www.altimed.co.uk)


Altiven compression hosiery is specifically designed for vascular conditions, lymphoedema and oedema. It is available in 7 different sizes and Made-to-Measure for a tailored fit. It is also available in a variety of styles:



**BELOW  
KNEE**



**THIGH  
LENGTH**




**THIGH WITH WAIST  
ATTACHMENT**



**HOLD UPS  
PLAIN**



**TIGHTS**



**HOLD UPS  
LACE**

### INDICATIONS

#### Low Compression (14 – 18 mmHg at the ankle)

Early onset and mild lymphoedema with minimal or no shape distortion. Stable condition. Elderly/arthritis patients. Pressure sensitive patients. ISL stages 0-II, maintenance, palliation, lipoedema, controlled cardiac oedema, dependency oedema, neurological deficit. May be considered for patients with ABPI 0.5-0.8 following full assessment. Prophylaxis.

#### European Class 1 (18 – 21 mmHg at the ankle)

Early onset and mild lymphoedema with minimal shape distortion. Lipoedema. Chronic Oedema. Stabilised condition. ISL stages 0-II, maintenance, palliation, lipoedema, controlled cardiac oedema, dependency oedema. ABPI >0.8.

#### European Class 2 (23 – 32 mmHg at the ankle)

Mild to Moderate lymphoedema, with minimal shape distortion. Maintenance after a leg ulcer has healed. Lipoedema. Long term maintenance. Varicose Veins. Venous conditions. ISL late stage II-III, phlebolympheoedema.

#### European Class 3 (34 – 46 mmHg at the ankle)

Moderate lymphoedema, with no shape distortion. To aid ulcer healing. Reduction therapy. Varicose Veins. Venous conditions. ISL stage III, phlebolympheoedema.



### CONTRAINDICATIONS


- Peripheral Vascular Disease (PVD), confirmed with an ABPI (Doppler) reading of less than 0.8
- Acute or unstable Heart Failure
- Ischaemic Heart Disease (IHD), also called Coronary Artery Disease
- Untreated Cellulitis or Phlebitis of the Lower Limb

#### Proceed with caution in those patients who have:

- Diabetes • Rheumatoid Arthritis

These patients may require referral to a specialist for advice and additional assessment. Close supervision may be necessary in line with local policies/pathways.

MEASURING AND FITTING GUIDE




Measurements in cms

Size	I	II	III	IV	V	VI	VII
G Hold-up	43-48	45-52	49-56	53-60	56-64	60-68	64-72
G Hold-up Ex Wide	49-57	53-62	57-67	61-72	65-77	69-81	73-85
G	43-57	45-62	49-67	53-72	56-77	60-81	64-85
D	27-33	29-36	32-39	34-42	36-45	38-48	40-50
C	28-34	30-37	33-40	35-43	37-46	39-49	41-51
B	18-20	20-22	22-24	24-26	26-28	28-30	30-32
Thigh (LG)				Below Knee (LD)			
Regular			72-83	Regular		39-44	
Short			62-71	Short		34-38	

EUROPEAN STANDARD ARMSLEEVES

Altiven armsleeves are designed for the specific requirements of lymphoedema. Improved stretch makes the sleeve easier to put on. They are made specifically to each patient's measurements and provide optimum compression and graduation along the entire arm. Altiven offers exceptional comfort with a smooth, soft and lightweight feel.



INDICATIONS

Low Compression

(12 - 16 mmHg)

Early onset and mild lymphoedema with minimal or no shape distortion. Stable condition. Elderly/arthritic patients. Pressure sensitive patients. Night maintenance. ISL stage I-II Prevention.



Class 1

(18 - 21 mmHg)

Mild lymphoedema with minimal shape distortion. Stable condition. Long term maintenance. ISL stage I-II.

Class 2

(23 - 32 mmHg)

Mild/Moderate lymphoedema, with minimal shape distortion. Long term maintenance. ISL stage II.



To order our Made-to-Measure hosiery, please request an order pad or order online at [www.altimed.co.uk](http://www.altimed.co.uk)

**Sockaid™**

Developed by patients, Sockaid is designed to simplify sock, stocking and compression hosiery fitting.

BENEFITS OF A SOCKAID:

- ✓ Easy to use and easy to follow instructions
- ✓ Can be used in or out of bed
- ✓ Lightweight and robust

Available on FP10/GP10



# Altipress®

## ALTIPRESS 40 HOSIERY KITS

Kits comprise 1 Medical Stocking & 2 liners

### SHORT LENGTH

MPC CODE	PIP CODE	SIZE	COLOUR
55SKO12	337-1325	S	●
55SKO22	337-1333	M	●
55SKO32	337-1341	L	●
55SKO42	337-1358	XL	●
55SKO52	337-1366	XXL	●

### REGULAR LENGTH

55RKO12	337-1374	S	●
55RKO22	337-1382	M	●
55RKO32	337-1390	L	●
55RKO42	337-1408	XL	●
55RKO52	337-1416	XXL	●

### LONG LENGTH

55LKO12	337-1424	S	●
55LKO22	337-1432	M	●
55LKO32	337-1440	L	●
55LKO42	337-1457	XL	●
55LKO52	337-1465	XXL	●

## ALTIPRESS LINER PACKS

Liner Packs contain 3 liners

### SHORT LENGTH

MPC CODE	PIP CODE	SIZE	COLOUR
55SLC10	337-1473	S	○
55SLC20	337-1481	M	○
55SLC30	337-1499	L	○
55SLC40	337-1507	XL	○
55SLC50	337-1515	XXL	○
55SLC12	341-8506	S	●
55SLC22	341-8696	M	●
55SLC32	341-8688	L	●
55SLC42	341-8670	XL	●
55SLC52	341-8613	XXL	●

### REGULAR LENGTH

55RLC10	337-1523	S	○
55RLC20	337-1531	M	○
55RLC30	337-1549	L	○
55RLC40	337-1556	XL	○
55RLC50	337-1564	XXL	○
55RLC12	341-8662	S	●
55RLC22	341-8654	M	●
55RLC32	341-8647	L	●
55RLC42	341-8621	XL	●
55RLC52	341-8605	XXL	●

### LONG LENGTH

55LLC10	337-1572	S	○
55LLC20	337-1580	M	○
55LLC30	337-1598	L	○
55LLC40	337-1606	XL	○
55LLC50	337-1614	XXL	○
55LLC12	341-8589	S	●
55LLC22	341-8563	M	●
55LLC32	341-8530	L	●
55LLC42	341-8522	XL	●
55LLC52	341-8514	XXL	●



To order our  
Made-to-Measure hosiery,  
request an order pad or order  
online at [www.altimed.co.uk](http://www.altimed.co.uk)

# Altiform®

## CLASS I. BELOW KNEE: Closed Toe

MPC CODE	PIP CODE	SIZE	COLOUR
511KC11	296-2611	S	●
511KC21	296-2629	M	●
511KC31	296-2637	L	●
511KC41	296-2645	XL	●
511KC12	296-2736	S	●
511KC22	296-2744	M	●
511KC32	296-2751	L	●
511KC42	296-2769	XL	●

## CLASS I. BELOW KNEE: Open Toe

511KO12	296-2694	S	●
511KO22	296-2702	M	●
511KO32	296-2710	L	●
511KO42	296-2710	XL	●

## CLASS I. THIGH LENGTH: Closed Toe

511TC11	296-2652	S	●
511KT21	296-2660	M	●
511KT31	296-2678	L	●
511KT41	296-2686	XL	●
511KT12	296-2819	S	●
511KT22	296-2827	M	●
511KT32	296-2835	L	●
511KT42	296-2843	XL	●

## CLASS I. THIGH LENGTH: Open Toe

511TO12	296-2777	S	●
511TO22	296-2785	M	●
511TO32	296-2793	L	●
511TO42	296-2801	XL	●

## CLASS II. BELOW KNEE: Closed Toe

MPC CODE	PIP CODE	SIZE	COLOUR
512KC11	296-2892	S	●
512KC21	296-2900	M	●
512KC31	296-2918	L	●
512KC41	296-2926	XL	●
512KC12	296-3015	S	●
512KC22	296-3023	M	●
512KC32	296-3031	L	●
512KC42	296-3049	XL	●

## CLASS II. BELOW KNEE: Open Toe

512KO12	296-2975	S	●
512KO22	296-2983	M	●
512KO32	296-2991	L	●
512KO42	296-3007	XL	●

## CLASS II. THIGH LENGTH: Closed Toe

512TC11	296-2934	S	●
512TC21	296-2942	M	●
512TC31	296-2959	L	●
512TC41	296-2967	XL	●
512TC12	296-3098	S	●
512TC22	296-3106	M	●
512TC32	296-3114	L	●
512TC42	296-3122	XL	●

## CLASS II. THIGH LENGTH: Open Toe

512TO12	296-3056	S	●
512TO22	296-3064	M	●
512TO32	296-3072	L	●
512TO42	296-3080	XL	●

## CLASS III. BELOW KNEE: Closed Toe

MPC CODE	PIP CODE	SIZE	COLOUR
513KC12	296-3338	S	●
513KC22	296-3346	M	●
513KC32	296-3353	L	●
513KC42	296-3361	XL	●

## CLASS III. BELOW KNEE: Open Toe

513KO12	296-3296	S	●
513KO22	296-3304	M	●
513KO32	296-3312	L	●
513KO42	296-3320	XL	●

## CLASS III. THIGH LENGTH: Closed Toe

513TC12	296-3411	S	●
513TC22	296-3429	M	●
513TC32	296-3437	L	●
513TC42	296-3445	XL	●

## CLASS III. THIGH LENGTH: Open Toe

513TO12	296-3379	S	●
513TO22	296-3387	M	●
513TO32	296-3395	L	●
513TO42	296-3403	XL	●



To order our  
Made-to-Measure hosiery,  
request an order pad or order  
online at [www.altimed.co.uk](http://www.altimed.co.uk)

## CLASS 1. BELOW KNEE: Size I

CODE	TOE	LENGTH	COLOUR
541K1C1R	Closed	Regular	●
541K1C2R	Closed	Regular	●
541K1C1S	Closed	Short	●
541K1C2S	Closed	Short	●
541K1O1R	Open	Regular	●
541K1O2R	Open	Regular	●
541K1O1S	Open	Short	●
541K1O2S	Open	Short	●

## CLASS 1. BELOW KNEE: Size II

541K2C1R	Closed	Regular	●
541K2C2R	Closed	Regular	●
541K2C1S	Closed	Short	●
541K2C2S	Closed	Short	●
541K2O1R	Open	Regular	●
541K2O2R	Open	Regular	●
541K2O1S	Open	Short	●
541K2O2S	Open	Short	●

## CLASS 1. BELOW KNEE: Size III

541K3C1R	Closed	Regular	●
541K3C2R	Closed	Regular	●
541K3C1S	Closed	Short	●
541K3C2S	Closed	Short	●
541K3O1R	Open	Regular	●
541K3O2R	Open	Regular	●
541K3O1S	Open	Short	●
541K3O2S	Open	Short	●

## CLASS 1. BELOW KNEE: Size IV

541K4C1R	Closed	Regular	●
541K4C2R	Closed	Regular	●
541K4C1S	Closed	Short	●
541K4C2S	Closed	Short	●
541K4O1R	Open	Regular	●
541K4O2R	Open	Regular	●
541K4O1S	Open	Short	●
541K4O2S	Open	Short	●

## CLASS 1. BELOW KNEE: Size V

541K5C1R	Closed	Regular	●
541K5C2R	Closed	Regular	●
541K5C1S	Closed	Short	●
541K5C2S	Closed	Short	●
541K5O1R	Open	Regular	●
541K5O2R	Open	Regular	●
541K5O1S	Open	Short	●
541K5O2S	Open	Short	●

## CLASS 1. BELOW KNEE: Size VI

541K6C1R	Closed	Regular	●
541K6C2R	Closed	Regular	●
541K6C1S	Closed	Short	●
541K6C2S	Closed	Short	●
541K6O1R	Open	Regular	●
541K6O2R	Open	Regular	●
541K6O1S	Open	Short	●
541K6O2S	Open	Short	●

## CLASS 1. BELOW KNEE: Size VII

541K7C1R	Closed	Regular	●
541K7C2R	Closed	Regular	●
541K7C1S	Closed	Short	●
541K7C2S	Closed	Short	●
541K7O1R	Open	Regular	●
541K7O2R	Open	Regular	●
541K7O1S	Open	Short	●
541K7O2S	Open	Short	●

## CLASS 2. BELOW KNEE: Size I

CODE	TOE	LENGTH	COLOUR
542K1C1R	Closed	Regular	●
542K1C2R	Closed	Regular	●
542K1C1S	Closed	Short	●
542K1C2S	Closed	Short	●
542K1O1R	Open	Regular	●
542K1O2R	Open	Regular	●
542K1O1S	Open	Short	●
542K1O2S	Open	Short	●

## CLASS 2. BELOW KNEE: Size II

542K2C1R	Closed	Regular	●
542K2C2R	Closed	Regular	●
542K2C1S	Closed	Short	●
542K2C2S	Closed	Short	●
542K2O1R	Open	Regular	●
542K2O2R	Open	Regular	●
542K2O1S	Open	Short	●
542K2O2S	Open	Short	●

## CLASS 2. BELOW KNEE: Size III

542K3C1R	Closed	Regular	●
542K3C2R	Closed	Regular	●
542K3C1S	Closed	Short	●
542K3C2S	Closed	Short	●
542K3O1R	Open	Regular	●
542K3O2R	Open	Regular	●
542K3O1S	Open	Short	●
542K3O2S	Open	Short	●

## CLASS 2. BELOW KNEE: Size IV

542K4C1R	Closed	Regular	●
542K4C2R	Closed	Regular	●
542K4C1S	Closed	Short	●
542K4C2S	Closed	Short	●
542K4O1R	Open	Regular	●
542K4O2R	Open	Regular	●
542K4O1S	Open	Short	●
542K4O2S	Open	Short	●

## CLASS 2. BELOW KNEE: Size V

542K5C1R	Closed	Regular	●
542K5C2R	Closed	Regular	●
542K5C1S	Closed	Short	●
542K5C2S	Closed	Short	●
542K5O1R	Open	Regular	●
542K5O2R	Open	Regular	●
542K5O1S	Open	Short	●
542K5O2S	Open	Short	●

## CLASS 2. BELOW KNEE: Size VI

542K6C1R	Closed	Regular	●
542K6C2R	Closed	Regular	●
542K6C1S	Closed	Short	●
542K6C2S	Closed	Short	●
542K6O1R	Open	Regular	●
542K6O2R	Open	Regular	●
542K6O1S	Open	Short	●
542K6O2S	Open	Short	●

## CLASS 2. BELOW KNEE: Size VII

542K7C1R	Closed	Regular	●
542K7C2R	Closed	Regular	●
542K7C1S	Closed	Short	●
542K7C2S	Closed	Short	●
542K7O1R	Open	Regular	●
542K7O2R	Open	Regular	●
542K7O1S	Open	Short	●
542K7O2S	Open	Short	●

## CLASS 1. STANDARD HOLD-UPS: Size I

CODE	TOE	LENGTH	COLOUR
541E1C1R	Closed	Regular	●
541E1C2R	Closed	Regular	●
541E1O1R	Open	Regular	●
541E1O2R	Open	Regular	●
541E1C1S	Closed	Short	●
541E1C2S	Closed	Short	●
541E1O1S	Open	Short	●
541E1O2S	Open	Short	●

## CLASS 1. STANDARD HOLD-UPS: Size II

541E2C1R	Closed	Regular	●
541E2C2R	Closed	Regular	●
541E2O1R	Open	Regular	●
541E2O2R	Open	Regular	●
541E2C1S	Closed	Short	●
541E2C2S	Closed	Short	●
541E2O1S	Open	Short	●
541E2O2S	Open	Short	●

## CLASS 1. STANDARD HOLD-UPS: Size III

541E3C1R	Closed	Regular	●
541E3C2R	Closed	Regular	●
541E3O1R	Open	Regular	●
541E3O2R	Open	Regular	●
541E3C1S	Closed	Short	●
541E3C2S	Closed	Short	●
541E3O1S	Open	Short	●
541E3O2S	Open	Short	●

## CLASS 1. STANDARD HOLD-UPS: Size IV

541E4C1R	Closed	Regular	●
541E4C2R	Closed	Regular	●
541E4O1R	Open	Regular	●
541E4O2R	Open	Regular	●
541E4C1S	Closed	Short	●
541E4C2S	Closed	Short	●
541E4O1S	Open	Short	●
541E4O2S	Open	Short	●

## CLASS 1. STANDARD HOLD-UPS: Size V

541E5C1R	Closed	Regular	●
541E5C2R	Closed	Regular	●
541E5O1R	Open	Regular	●
541E5O2R	Open	Regular	●
541E5C1S	Closed	Short	●
541E5C2S	Closed	Short	●
541E5O1S	Open	Short	●
541E5O2S	Open	Short	●

## CLASS 1. STANDARD HOLD-UPS: Size VI

541E6C1R	Closed	Regular	●
541E6C2R	Closed	Short	●
541E6O1R	Open	Regular	●
541E6O1S	Open	Short	●
541E6C2S	Closed	Regular	●
541E6C2S	Closed	Short	●
541E6O2R	Open	Regular	●
541E6O2S	Open	Short	●

## CLASS 1. STANDARD HOLD-UPS: Size VII

541E7C1R	Closed	Regular	●
541E7C1S	Closed	Short	●
541E7O1R	Open	Regular	●
541E7O1S	Open	Short	●
541E7C2R	Closed	Regular	●
541E7C2S	Closed	Short	●
541E7O2R	Open	Regular	●
541E7O2S	Open	Short	●

### CLASS 2. STANDARD HOLD-UPS: Size I

CODE	TOE	LENGTH	COLOUR
542E1C1R	Closed	Regular	●
542E1C1S	Closed	Short	●
542E1O1R	Open	Regular	●
542E1O1S	Open	Short	●
542E1C2R	Closed	Regular	●
542E1C2S	Closed	Short	●
542E1O2R	Open	Regular	●
542E1O2S	Open	Short	●

### CLASS 2. STANDARD HOLD-UPS: Size II

542E2C1R	Closed	Regular	●
542E2C1S	Closed	Short	●
542E2O1R	Open	Regular	●
542E2O1S	Open	Short	●
542E2C2R	Closed	Regular	●
542E2C2S	Closed	Short	●
542E2O2R	Open	Regular	●
542E2O2S	Open	Short	●

### CLASS 2. STANDARD HOLD-UPS: Size III

542E3C1R	Closed	Regular	●
542E3C1S	Closed	Short	●
542E3O1R	Open	Regular	●
542E3O1S	Open	Short	●
542E3C2R	Closed	Regular	●
542E3C2S	Closed	Short	●
542E3O2R	Open	Regular	●
542E3O2S	Open	Short	●

### CLASS 2. STANDARD HOLD-UPS: Size IV

542E4C1R	Closed	Regular	●
542E4C1S	Closed	Short	●
542E4O1R	Open	Regular	●
542E4O1S	Open	Short	●
542E4C2R	Closed	Regular	●
542E4C2S	Closed	Short	●
542E4O2R	Open	Regular	●
542E4O2S	Open	Short	●

### CLASS 2. STANDARD HOLD-UPS: Size V

542E5C1R	Closed	Regular	●
542E5C1S	Closed	Short	●
542E5O1R	Open	Regular	●
542E5O1S	Open	Short	●
542E5C2R	Closed	Regular	●
542E5C2S	Closed	Short	●
542E5O2R	Open	Regular	●
542E5O2S	Open	Short	●

### CLASS 2. STANDARD HOLD-UPS: Size VI

542E6C1R	Closed	Regular	●
542E6C1S	Closed	Short	●
542E6O1R	Open	Regular	●
542E6O1S	Open	Short	●
542E6C2R	Closed	Regular	●
542E6C2S	Closed	Short	●
542E6O2R	Open	Regular	●
542E6O2S	Open	Short	●

### CLASS 2. STANDARD HOLD-UPS: Size VII

542E7C1R	Closed	Regular	●
542E7C1S	Closed	Short	●
542E7O1R	Open	Regular	●
542E7O1S	Open	Short	●
542E7C2R	Closed	Regular	●
542E7C2S	Closed	Short	●
542E7O2R	Open	Regular	●
542E7O2S	Open	Short	●

 To order our  
Made-to-Measure hosiery,  
request an order pad or order  
online at [www.altimed.co.uk](http://www.altimed.co.uk)

**Sockaid™**

MPC CODE	PIP CODE
60000	319-4016

To find out more please call us on:  
**01509 501720**  
or visit us online at  
[www.urgo.co.uk](http://www.urgo.co.uk) [www.altimed.co.uk](http://www.altimed.co.uk)

Altipress® Altiform® Altiven®

**URGO**  
MEDICAL  
Healing people®



MD 648377



Accredits Urgo Medical research in  
wound healing and healthy skin

Urgo Limited - Sullington Road, Shepshed, Loughborough, LE12 9JG

\*Manufactured and dispatched within 72 hours of receipt of the order. Notwithstanding forces beyond our control.  
Altiform, Altipress and Altiven are all registered trademarks of Urgo Limited

**MADE IN  
BRITAIN**